



Cross-Linked Hyaluronic Acid Filler

For Internal training

contents

01. Characteristics of Revitrane

02. Revitrane Product Range

0.3 Before & After

04. Certificates and Facilities



Revitrane is..

Revitrane – Maximized advantage of HA filler

Safe monophasic type filler with minimal crosslinking agent(BDDE).

Safety

- High-quality raw material
- Minimized chemical cross-linking rate using BDDE
- Advanced residual BDDE cleaning process

Excellent Viscoelasticity

- Uniform cross-linked HA structure
- Soft viscoelastic gel texture for making natural volume

Cohesivity

- Keep initial volume for long time by cohesive gel
- Improved skin integration that blends with surrounding skin tissue

Soft Injection

- Smooth and uniform injection
- Low extrusion force compare to the others

Revitrane is..

Monophasic type

Soft viscoelastic hyaluronic acid gel texture



REVITRANE

- *Classic*
- *Volume*
- *Premium*

CE
2292

Revitrane is..

Skin Booster


Cross-linked gel with smooth and low cohesion gel for biorevitalization and moisturizing

REVITRANE

HA20 Skin Booster



Safety



Hyaluronic Acid

High Quality HA Raw materials

Approved raw material (from Shiseido, Japan)

Revitrane is manufactured using safe and qualified hyaluronic acid approved US FDA and European Directorate for the Quality of Medicines(EDQM)

Constant HA molecular structure

Stable HA structure that can endure external pressure

Purity


Safe material with less than 10ppm of heavy metal (standards for US FDA and EDQM is 20ppm)

	PH. Eur	Revitrane
Residual protein	Maximum 0.3% or 0.1%	Maximum 0.1%
Endotoxin	<0.5 or 0.05IU/mg	<0.04IU/mg

Safety

BR PHARM's own HA cross-linking process

Low-chemical & Enhanced Double Cross-linking Method : **LED Method**



Minimizes incomplete crosslinking parts(Pendant cross-linker) and increase complete cross-linking rate using long-time and low-temperature cross-linking process.

Safety


Purification Technology of Cross-Linked HA

Advanced Technique for Elimination of Residues, ATER

Unique washing technology removes residual BDDE and impurities with minimized HA deformation



Factor	J Product	Revitrane
BDDE amount	800~1000ppm	500ppm



Color comparison after drying HA gel




* MOD : Degree of modification

Viscoelasticity

Excellent Viscoelasticity(G')

Revitrane is soft gel type filler with very fine HA particles for optimal viscoelasticity makes natural volume. And Revitrane is easily molded and smooth injection feeling.

Product	0.1Hz	1.0Hz
Juvederm Voluma	310~330 Pa	350~400 Pa
Revitrane Premium	328 Pa	383 Pa



Revitrane Premium has similar viscoelasticity to Juvederm Voluma product

Cohesivity

High Cohesivity

Revitrane's cohesive gel texture maintains initial shape and minimizes migration under pressure

Ink diffusion testing result after 6 hours for Revitrane Premium and other Korean products



**Revitrane
Premium**



A

B


Soft Injection

Low extrusion force improves treatment satisfaction

Revitrane has low and constant extrusion force though high-quality HA material and purification process for making uniform HA structure.

Product	Extrusion force(N)
Juvederm Voluma	12~15
Revitrane Premium	10~12

Extrusion Force vs Distance




Soft Injection


Low extrusion force improves treatment satisfaction

Revitrane Premium, Volume, Classic has an low injection pressure of 12 to 19N when using a 30G needle and an injection pressure of 3 to 7N when using a 27G needle. It is possible to select the needle according to the skill of the operator with a lower injection pressure than other filler products to perform the procedure.


Revitrane Classic



Revitrane Volume



Revitrane Premium



Revitrane – Biphasic type lines

Revitrane Premium Q and Premium SUB-Q

Particle structured gel type with high elasticity for extensive facial volume augmentation



Revitrane – Biphasic type lines

Biphasic type

Uniform cross-linked hyaluronic acid particles

- **Stabilizing & Enhancing Elasticity Double-cross-linking(SEED) Technology**
- **Maximize binding ratio by our own cross-linking method for biphasic structure**
 - Chemical crosslinks and twist structure
 - Control and make different particle size for treating purpose

- **Control the size of particles**

Product	Particle Size(micron)
Revitrane Premium SUB-Q	800
Revitrane Premium-Q	600

Particles (SUBQ) ►



contents

01. Characteristics of Revitrane

02. Revitrane Product Range









03. Before & After

04. Certificates and Facilities








Product Specification & Treatment Area

Revitrane Series (Monophasic)

Product	Packaging	Concentration	Needle	Injection depth
 Skinbooster	2ml*3syr/box	20mg/ml	31G or 34G(NanoNeedle) Multi Needle Injector	
 Classic	1ml*1syr/box	20mg/ml	30G	
 Volume	1ml*1syr/box	20mg/ml	27G	
 Premium	1ml*1syr/box	20mg/ml	27G	

Product Specification & Treatment Area

Revitrane Premium Series(Biphasic)

Product	Packaging	Concentration	Needle	Injection depth
 Revitrane Premium Q	1ml*1syr/box	20mg/ml	25G	
 Revitrane Premium SUB-Q				
 Revitrane Premium SUB-Q	1ml*1syr/box	20mg/ml	25G	

Revitrane HA20 Skin Booster

Revitrane HA20 Skin Booster

Skin Booster has lower viscosity and elasticity, and is suitable for improving fine wrinkles and overall moisturizing & elasticity of the face.



Gel Type

- Viscosity ●
- Elasticity ●
- Cohesivity ●

Needle Size

- Needle is not included
- 31G is suitable


Treatment Area

●
Skin
Booster

- Whole face
- Water shine treatment



Injection Depth



Revitrane *Classic*

Revitrane Classic

Classic is designed to be treatment of the areas like lips or fine lines for expressing natural volume. It is suitable for areas that require natural volume due to less swelling.

Gel Type

- Viscosity ● ●
- Elasticity ●
- Cohesivity ●

Treatment Area

●
Classic


- Around eyes
- Under eyes
- Smile lines
- Wrinkles around lips
- Lip volume



Needle Size

- 30G x 2 needles are included
- 30G-31G needle is suitable

Injection Depth



Revitrane Volume

Revitrane Volume

Volume is viscous gel texture to restore volume as it attracts and keeps moisture in the skin.

Volume has high expansion rate, so it is suitable for volumizing

Gel Type

- Viscosity ● ● ●
- Elasticity ●
- Cohesivity ● ●

Treatment Area

●
Volume


- Cheek
- Nasolabial folds



Needle Size

- 27G x 2 needles are included
- 27-30G needle is suitable

Injection Depth



Revitrane Premium

Revitrane Premium

Premium is slightly viscous and high elastic gel type, so it is suitable for most of treatment areas for volumizing



Gel Type

- Viscosity ●
- Elasticity ● ● ●
- Cohesivity ● ● ●


Needle Size

- 27G x 2 needles are included
- 27-30G needle is suitable


Treatment Area

●
Premium

- Forehead
- Cheek bone
- Nose
- Nasolabial folds
- Chin



Injection Depth



Revitrane *Premium-Q*

Revitrane Premium-Q

Premium-Q's biphasic type particles provides good support to maintain structure, it is suitable for making shape of the nose area.



Gel Type

- Viscosity -
- Elasticity ● ● ●
- Cohesivity ●


Needle Size

- 25G x 2 needles are included


Treatment Area

●
Premium

- Nose tip
- Dorsum of the nose



Injection Depth



Revitrane Premium SUB-Q

Revitrane Premium SUB-Q

Premium SUB-Q has bigger size biphasic type particles than Premium-Q and it provides good support to maintain structure for cheek or chin area



Gel Type

- Viscosity -
- Elasticity ● ● ●
- Cohesivity ●

Needle Size


- 25G x 2 needles are included

Treatment Area

- Premium SUB-Q
- Cheek augmentation
- Chin augmentation



Injection Depth



contents

01. Characteristics of Revitrane

02. Revitrane Product Range

0.3 Before & After

04. Certificates and Facilities



Clinical data

Before & After

Revitrane Volume, Age 43, Female, Nasolabial fold Left 2cc/Right 3cc

Before



Right After



After (7days later)



Clinical data

Before & After

Revitrane Volume, Age 40, Female, Nasolabial fold Left 3cc/Right 3cc

Before



Right After



After (7days later)



Clinical data

Before & After

Lips : Revitrane Classic 0.8ml



Clinical data

Before & After



1. Cheek : **Revitrane Premium SUB-Q 1cc**
2. Whole face : **Revitrane HA20 Skin Booster**



Clinical data

Before & After

1. Forehead + eyebrow
: **Revitrane Classic 1.2ml**
2. 1st and 2nd zygomatic lines
: **Revitrane Volume 1.0ml**
3. Nasolacrimal groove
: **Revitrane Classic 0.6ml**
4. Nasolabial folds
: **Revitrane Classic 1.4ml**
5. Lips
: **Revitrane Classic 0.8ml**
6. Whole face
: **Revitrane Skin booster 2.0ml**



contents

01. Characteristics of Revitrane

02. Revitrane Product Range

03. Before & After

04. Certificates and Facilities



Study

Pre-clinical study



In the comparative advantage test, as a result of testing the volume of Revitrane and J Company's product in the folds, the Rebitran volume showed excellent initial volume and durability.

In a comparative study on injections between Revitrane volume and Product J, The result showed great efficacy and long lasting at the nasolabial folds

- **Efficacy assessment – Long lasting Scored by Chung Ang Univ., Korea 2018**

Volume Change after Revitrane

Volume injection



The change of initial filler volume: Volumetric analysis; PRIMOS_(LITE) topography (with software PRIMOS 5.8)

Study

Research of cross-linking #1

Study of swelling degree and elasticity according to HA molecular weight

Polymer(Korea), Vol. 39, No. 6, pp. 976-980 (2015)
<http://dx.doi.org/10.7317/pk.2015.39.6.976>

ISSN 0379-153X(Print)
 ISSN 2234-8077(Online)

히알루론산 피부용 필러의 분자량에 따른 팽윤도와 탄성계수 영향

이득용^{*} · 전철병^{*} · 손시원^{*} · 김영주^{**} · 김진태^{**} · 장주용^{***} · 김석순^{****}

대림대학교 의공융합과, ^{*}베리콤 기술연구소, ^{**}네오바이오텍 연구소,

^{***}셀부메드 생체재료의공학연구소 ^{****}유비알뷰티플 레볼루션

(2015년 7월 16일 접수, 2015년 7월 23일 수정, 2015년 8월 3일 채택)

Influence of Molecular Weight on Swelling and Elastic Modulus of Hyaluronic Acid Dermal Fillers

Deuk Yong Lee¹, Cheolbyung Cheon, Siwon Son, Young-Zu Kim^{2*}, Jin-Tae Kim^{2*},
 Ju-Woong Jang^{3**} and Seok-Soon Kim^{4**}

¹Department of Biomedical Engineering, Daelim University, Anyang 431-715, Korea

²R&D Center, Vericom Co., Ltd., Anyang 606-72, Korea

³R&D Center, Neobiotech Co., Ltd., Seong 152-789, Korea

⁴R&D Laboratory, Cellumed Co., Ltd., Seoul 153-782, Korea

⁵Beautiful Revolution Co., Ltd., Seoul 135-513, Korea

(Received July 16, 2015; Revised July 23, 2015; Accepted August 3, 2015)

초록: 분자량이 다른 생분해성 히알루론산(HA)에 1,4-butanediol diglycidyl ether(BDDE)가 가교된 HA를 첨가하여 히드로겔을 제조하고 HA 분자량에 따른 팽윤도와 탄성계수의 효과를 조사하였다. 최종 무석 후 잔류 BDDE 양은 0.5 ppm 이하이었다. 비가교된 분자량이 큰 HA에서 최대 팽윤도가 관찰되었다. 팽윤도는 비가교 HA 양과 가교도에 각각 비례와 반비례하였다. 1.0 w/v% BDDE 가교 HA와 HA(1368 kDa)를 혼합한 monophasic 히드로겔의 탄성계수는 152~325 Pa이었다. 15% 가교 HA(687 kDa)과 85% HA(1368 kDa)로 구성된 biphasic 히드로겔의 탄성계수는 178 Pa이었다.

978 Deuk Yong Lee, et al.

Results and Discussion

These different HA were crosslinked with BDDE (0.1, 0.5, 1.0 w/v%) to extend the distance between crosslinked HA to completely immobilized a few day after operation. The epoxide groups on BDDE preferentially react with the alcohol on the HA backbone forming an ether bond connection, as depicted in Figure 1. The presence of residual BDDE after addition was evaluated by using a GC, as shown in Figure 2. Since the BDDE, detected at 4.8 min in GC graph (Figure 2), is likely to be the unreacted product, the amount of unreacted

residual BDDE is maintained at trace amount of ~2 ppm, which has been determined to be safe after a Food and Drug Administration (FDA)'s GC results of 0.1 w/v% BDDE crosslinked HA. It is revealed that the amount of residual unreacted BDDE after final curing was less than 0.7 ppm, suggesting that the crosslinking was successfully occurred. The highest swelling ratio (SOR) is observed for the crosslinked HA (1:1) having the largest molecular weight (1368kDa), as depicted in Figure 3. The swelling ratio decreased with increasing the BDDE concentration and decreasing the molecular weight due to an increased number of rigid HA chain interaction. SEM results revealed that pore size and porosity decreased and increased with increasing the BDDE concentration, respectively, as illustrated in Figure 3. A 1/2 amount of HA(1) is formed when crosslinking between




Figure 1. Chemical structure of crosslinked HA with BDDE.




Figure 2. GC graph of 0.1 w/v% BDDE crosslinked HA (1:1).




Figure 3. Swelling ratio of various HA crosslinkers as a function of BDDE content.




Figure 4. SEM images of HA crosslinkers at HA (1:1), HA (1:2), and HA (1:3) crosslinkers.

"Epoxy resin" due to the reactive nature of the epoxide groups. FDA-approved HA is crosslinked with a design specification of ~2 ppm unreacted BDDE on the final product. GC results revealed that the amount of residual unreacted BDDE in HA(1) was maintained at trace amount of ~2 ppm, demonstrating that the BDDE-crosslinked HA(1) was clinically safe. Also, there remaining trace amount, which the FDA has determined to be below the level that is safe after a safety risk assessment, are easily converted into CO₂ and water. Since the filling capacity of a dermal filler is known to be dependent on its water uptake capacity, the swelling ratio was

significantly higher for monophasic crosslinked gel prepared by mixing the HA and crosslinking in one step, allowing for ease of injection. The presence of crosslinked HA in HA(1) are very suitable to Veribest products (80% crosslinker used in the market).¹⁴⁻¹⁶ However, Veribest is a biphasic product consisting of the crosslinked particles (77%) suspended in non-crosslinked HA (23%) used as a carrier. Staoplastic made from 87 kDa HA (1:1) was blended with 1368 kDa HA (1:1) in different ratio, as depicted in Figure 7. Heterophasic was crosslinked after the crosslinking of hydroxyl groups by this crosslinker.

Research of cross-linking #2

Polymer(Korea), Vol. 40, No. 4, pp. 600-606 (2016)
<http://dx.doi.org/10.7317/pk.2016.40.4.600>

ISSN 0379-153X(Print)
 ISSN 2234-8077(Online)

가교된 히알루론산 구슬로 제조한 피부용 필러의 점탄성 특성

전철병 · 김예나 · 손시원 · 이득용[†] · 김진태 · 권미경* · 김영주** · 김석순***
 대림대학교 의공융합과, *네오바이오텍 치과재료연구소, **베리콤 기술연구원, ***주비알뷰티플 레플루션
 (2016년 2월 2일 접수, 2016년 2월 26일 수정, 2016년 3월 7일 채택)

Viscoelasticity of Hyaluronic Acid Dermal Fillers Prepared by Crosslinked HA Microspheres

Cheolbyong Chun, Yena Kim, Siwon Son, Deuk Yong Lee[†], Jin-Tae Kim, Mi-Kyung Kwon*,
 Young-Zu Kim**, and Seok-Soon Kim***

Department of Biomedical Engineering, Daelim University, Anyang 13916, Korea
[†]R&D Center, Neobiotech Co., Ltd, Seoul 06881, Korea
 **R&D Center, Vericom Co., Ltd, Anyang 14087, Korea
 ***Beautiful Revolution Co., Ltd., Seoul 06151, Korea

(Received February 2, 2016; Revised February 26, 2016; Accepted March 7, 2016)

초록: 다비닐 실론 가교제로 가교한 히알루론산(HA) 구슬과 비가교된 HA의 부피비가 65/35-95/5로 다른 피부용 HA 필러를 제조하여 가교된 HA 구슬이 필러의 탄성계수와 입자감에 대한 효과를 조사하였다. 구슬 내 2-4±0.5 μm 내부기공을 가진 HA 구슬의 평균 입도는 60~100±4 μm이었다. HA 필러의 가교된 젤 입자 크기는 300±30 μm이었다. 가교된 HA 미세구슬의 부피비가 65에서 95%로 증가함에 따라 필러의 탄성계수는 211에서 700 Pa로 증가하였다. 29~30 게이지 주사바늘을 통과할 수 있는 175~420 Pa 탄성계수를 가진 필러의 가교된 구슬의 부피함량은 65~85%이었다. 실험결과, 모든 필러들은 가교된 HA 구슬의 부피비가 증가할수록 젤 입자 밀도 증가로 인하여 입자감도 증가하였다. 본 연구에서 injectability와 입자감이 우수한 피부용 필러를 성공적으로 제조하였다.

602

C. Chun et al.

preparation, free HA (20 mg/mL) dissolved in a phosphate buffered saline solution (PBS, NaH₂PO₄) and swollen cross-linked microspheres were homogenized for 1-4 min and then incubated for 24 h.⁶ The ratio of the crosslinked HA to the free HA is varied from 65:35 to 95:5. The elastic and viscous response of hydrogel depend on the concentration and molecular weight of the HA and on the frequency used during the measurements. Rheological behavior of HAHs were analyzed with a Thermo Haake RS1 Rheometer (Newington, USA), using a plate and plate geometry with a 1.2 mm gap.¹² All measurements were performed using a 20 mm titanium sensor at 25 °C. Oscillation measurements were taken at 5 Pa tau over a frequency range of 0.01 to 100 Hz. The measured properties of hydrogels was evaluated by measuring the elastic modulus (G') and viscous modulus (G''). The measured modulus at frequency of 5 Hz for HAHs was compared.¹³ In addition, PTF is achieved by observing or touching HAHs using an optical microscope (Pro Carescope, Somatech, Korea) or fingers, respectively.

Swelling Property. The swelling characteristics were measured by immersing weighed samples of dry microspheres for 24 h in PBS. The gels were sectioned through a 3.0 μm membrane filter (Advantec, Japan). The excess surface water in the

swollen gel was removed by blotting and then the swollen gel was weighed. After measuring the weight of the gels, the swelling ratio (S) was determined by using eq. (1).

$$S(\%) = \frac{W_s - W_d}{W_d} \times 100 \quad (1)$$

Where, W_s and W_d are the weight of the swollen gel and the dry gel, respectively.

Gas Chromatography (GC). DVS reacts with the primary alcohol groups in the HA backbone. The crosslinked microspheres were cleaned in water and ethanol to eliminate an unreacted residual crosslinker. The presence of DVS after cleaning may cause adverse, allergic reactions and potential neurotoxicity of the dermal fillers because they are used within the dermis for several months. The presence of the unreacted residual crosslinker in HAHs is evaluated by using gas chromatography (YL6100 GC, Younglin Co., Ltd, Yongin, Korea).¹⁴

Results and Discussion

HA microspheres were fabricated by using a modified spray method.^{12,17} The experimental procedure was described else-




Figure 1. SEM image of (a),(b) HA microspheres crosslinked for 24 h. Note that (c)-(d) images exhibit porous network inside the HA microspheres.

원리, 제40권 제4호, 2016년

A between 65:35 of 175 to 420 Pa, after allowing for the treatment of r and softer skin kin deformation.⁶ sphere is in the styrene product act.¹⁴ The elastic (Pa) were even Pa).¹² acial Restylene idk. The biphasic it of crosslinked JA. Our previous files of the dermally, which was contain a less-





Figure 3. (a) Elastic modulus as a function of angular frequency at 25 °C; (b) the composition data of dermal fillers. Note that the grid of the dotted line in (b) represents the easy delivery of the fillers through a 29-30 gauge needle.

원리, 제40권 제4호, 2016년

ner chain interactions and thereby lowering the viscosity. Viscous modulus (G'') is also called loss modulus because it describes the energy that is lost as viscous dissipation. The G'' value is a measure of the flow (rheological) properties for the fillers. Although G'' of 70 Pa was observed for the fillers consisting of 65:35, as depicted in Figure 4, the rest of fillers showed higher G'' values higher than 128 Pa. The G'' values of the as-prepared fillers (128 to 166 Pa) are similar to those of Restylene (119 Pa) and Perlane (125 Pa), which are categorized in the medium-viscosity and medium-elasticity group.¹⁸ No appreciable flowability is detected. In addition, PTF is suc-




Figure 4. Viscous modulus as a function of angular frequency at 25 °C.

원리, 제40권 제4호, 2016년

Certificates



✓ CE : Cross-Linked Hyaluronic acid Gel filler (Class III)
Expiry Date : 19.DEC.2022


✓ ISO 13485

✓ GMP

Test reports

Pre-clinical test reports

YOUR PARTNER FOR THE BEST QUALITY



TEST REPORT (Registered Copy)

81, Gyoysaek-ro, Gwancheon-si, Gyeonggi-do, 138-1, Rep. of KOREA TEL: 82-2-2164-0011 FAX: 82-2-2504-1038

Report No.: MKR-09022 Receipt Date: Jun. 5, 2015
 Client: Seo, Gyeong-soon Test Completion Date: Aug. 03, 2015
 BR Pharm Co. Ltd. LOT NO.: R191013P
 811 Medical Device Complex, 200, Gaeoposi-ro, Jijang-myeon, Wonju-si, Gangwon-do, Korea

Overall Judgement: Passed MFG. Date: 2015.06.01
 Sample: Graft/prosthesis, biomaterial EXP. Date: -

TEST RESULTS

TEST ITEM	UNIT	SAMPLE	CRITERIA	RESULT	ACCEPT. CR.	TEST METHOD
Injection-Systemity JEPT-2015-02-0013P		Refractive HA10	Pass. Guide: (Slightly) Passed.	Passed.	Yield: 100%	Technical Documents (See note)
Injection - Uppare hyaluronic acid identification		Refractive HA10	The refractive index of the sample is more than 1.375 (min. 1.375)	Confirmed.	Yield: 100%	Technical documents

Extraction Condition
 (LAL TEST) : LAL reagent water (3 mL / device/mL)
 (Ophthalmic test) : 4 g / 42 mL, 37 °C, 24 h. The extract was centrifuged (250 g, 10 min)

- Test period: Jun. 15, 2015 - Jul. 31, 2015

Annex #1: Picture of appearance or shape
 Attachment: Test Report

USAGE: For approval


NOTE: 1. The test results on this test report are only limited to the samples and sample names provided by the customer and KTR does not guarantee the quality of all products of the customer, and you can confirm the authenticity of the test report online (www.ktr.or.kr) or by using the QR code.
 2. This test report shall not be used for public relation, advertisement, lawsuit and any other purposes outside the scope of its defined usage.

Lee Do-gyun
Prepared by: Lee Do-gyun
E-mail: leedo@ktr.or.kr

Seung-Wooh Park
Reviewed by: Seung-Wooh Park
Technical Manager
Tel.: 1577-2281/2815 (D-40)


Aug. 04, 2015 Registered copy date: Jun. 12, 2015

Korea Testing & Research Institute
President *Choi Hyeongha*




QR Code to verify genuineness

1 of Total: 3 Pages



YOUR PARTNER FOR THE BEST QUALITY



TEST REPORT (Registered Copy)

81, Gyoysaek-ro, Gwancheon-si, Gyeonggi-do, 138-1, Rep. of KOREA TEL: 82-2-2164-0011 FAX: 82-2-2504-1038

Report No.: MKR-09022 Receipt Date: Jun. 5, 2015
 Client: Seo, Gyeong-soon Test Completion Date: Oct. 14, 2015
 BR Pharm Co. Ltd. LOT NO.: S21218P
 811 Medical Device Complex, 200, Gaeoposi-ro, Jijang-myeon, Wonju-si, Gangwon-do, Korea

Overall Judgement: Passed MFG. Date: 2015.06.01
 Sample: Graft/prosthesis, biomaterial EXP. Date: -

TEST RESULTS

TEST ITEM	UNIT	SAMPLE	CRITERIA	RESULT	ACCEPT. CR.	TEST METHOD
Impaction test (ISO 12 series, instrumented)		Refractive HA10	pass	(0) Passed	Yield: 100%	Technical Documents (See note)

(1) As a result of subcutaneous implantation test or test article in NUOJ rabbits for 12 weeks, the test article was calculated as 0.00 compared to control histological evaluation.
 Thus, the biocompatibility rating was evaluated as Non-irritant.

- Test period: Jun. 25, 2015 - Oct. 14, 2015

Annex #1: Picture of appearance or shape
 Attachment: Test report

USAGE: For approval


NOTE: 1. The test results on this test report are only limited to the samples and sample names provided by the customer and KTR does not guarantee the quality of all products of the customer, and you can confirm the authenticity of the test report online (www.ktr.or.kr) or by using the QR code.
 2. This test report shall not be used for public relation, advertisement, lawsuit and any other purposes outside the scope of its defined usage.

Lee Do-gyun
Prepared by: Lee Do-gyun
E-mail: leedo@ktr.or.kr

Lee Seung-young
Reviewed by: Lee Seung-young
Technical Manager
Tel.: 1577-2281/2815 (D-40)


Oct. 14, 2015 Registered copy date: Jun. 12, 2015

Korea Testing & Research Institute
President *Choi Hyeongha*



QR Code to verify genuineness

1 of Total: 1 Page



YOUR PARTNER FOR THE BEST QUALITY



TEST REPORT (Registered Copy)

81, Gyoysaek-ro, Gwancheon-si, Gyeonggi-do, 138-1, Rep. of KOREA TEL: 82-2-2164-0011 FAX: 82-2-2504-1038

Report No.: MKR-09024 Receipt Date: Jun. 10, 2015
 Client: Seo, Gyeong-soon Test Completion Date: Aug. 03, 2015
 BR Pharm Co. Ltd. LOT NO.: R191013P
 811 Medical Device Complex, 200, Gaeoposi-ro, Jijang-myeon, Wonju-si, Gangwon-do, Korea

Overall Judgement: Passed MFG. Date: 2015.06.01
 Sample: Graft/prosthesis, biomaterial EXP. Date: -

TEST RESULTS

TEST ITEM	UNIT	SAMPLE	CRITERIA	RESULT	ACCEPT. CR.	TEST METHOD
Injection-Systemity Dissolution/In vitro test		Refractive HA10	pass	(3) Passed	Yield: 100%	Technical documents
Injection-Axide Systemic Toxicity Test		Refractive HA10	pass	(3) Passed	Yield: 100%	Technical documents
Injection-Systemity test(OA test)		Refractive HA10	pass	(3) Passed	Yield: 100%	Technical documents
Injection-Systemity Test (Subacute Toxicity Test)		Refractive HA10	Pass	Negative	Yield: 100%	Technical Documents
Injection-Systemity Test (Subacute Toxicity Test)		Refractive HA10	Pass	Negative	Yield: 100%	Technical Documents

- Next Page -

Jeong Seung-uk
Prepared by: Jeong Seung-uk
E-mail: jeungs@ktr.or.kr

Lee Seung-young
Reviewed by: Lee Seung-young
Technical Manager
Tel.: 1577-2281/2815 (D-40)

Aug. 03, 2015 Registered copy date: Jun. 12, 2015

Korea Testing & Research Institute
President *Choi Hyeongha*



QR Code to verify genuineness

1 of Total: 2 Pages



Facilities

Automated manufacturing facilities

Compliant with GMP & ISO13485 production process

- ✓ Clean room condition
 - controlled by ISO14644-1
- ✓ 15,000 ~ 20,000 syr/day production capacity
- ✓ GMP & ISO13485



Thank You

A hand in a white lab coat is holding a syringe and a molecular model. The syringe is labeled 'BR PHARM' and '10ml'. The molecular model consists of white, black, blue, and red spheres connected by white rods. The background is a blurred image of a person in a white lab coat.

BR PHARM CO., LTD. www.brpharm.com
E-mail : brpharm3@brpharm.com

Head Office & Factory

13, Sinpyeong-ro, Sinpyeong-ri, jjeong-myeon,
Wonju-si, Gangwon-do, Republic of Korea
Tel : +82-33-901-4500 , Fax : +82-33-901-4600

Research Institute & Overseas Business Div.

802, 803, C , Hyundai Knowledge Industry Center,7,
Beobwon-ro 11-gil, Songpa-gu, Seoul, Republic of Korea
Tel : +82-2-2148-6770 , Fax : +82-2-2148-6777

BR PHARM