

Pure Eyes

PN Booster



Pure Eyes

Pure Eyes is a PN booster that enables you to regain your natural beauty through hydration and collagen regeneration of the skin.

Product Name	Pure Eyes
Package	1.1ml x 1 syringe
Storage	Keep away from direct sunlight. Keep at room temperature (1~30°C)
Expiration date	24 months from the date of manufacture
Composition	PN (Polynucleotide) 2%

What is PN ?

PN stands for Polynucleotide.

It is extracted from salmon and it accelerates the generation of collagen and elastin which increases the density and thickness of the skin.

Premium medical grade PN is used for Pure Eyes.




SALMON



PN



PURE EYES



Benefits of PN

Pure Eyes PN Booster



SKIN REPAIR

Boosts cell regeneration.
Stimulates collagen production.




REJUVENATION

Promotes a more youthful and rejuvenated complexion by boosting collagen levels.



HYDRATION

Works by attracting and retaining moisture in the skin, improving hydration.



ELASTICITY IMPROVEMENT

Stimulates elastin to increase skin elasticity.



SKIN TONE IMPROVEMENT


Improves pigmentation and skin tone by increasing the thickness and density of the dermis.



PORE TIGHTENING

Provides pore tightening effects by filling open pores from within the skin.

Pure Eyes Indication



Powerful
Hydration



Elasticity
Improvement



Balance
Recovery



Skin
Regeneration

Treatment Guide

Purpose	To improve dark circle
Needle Size	32G 4mm
Technique	Micro droplet technique
Treatment	Apply 5-10 drops under one eye, (depending on the size of the area of application)
Recommended Session	3 Sessions (2-3 weeks interval)

Before & After



Pure Eyes

PN Booster

"Pure"

: without any extraneous and unnecessary elements.

(Oxford Languages)

www.hjcorporations.kr

info@hjmediclub.com

+82 2 2641 0209

 HJ corporations